Etapa 4 D Instale o cabo de Ethernet.

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GE Power
Water & Process Technologies

IMPORTANT NOTICE

For more detailed instructions, you can download manuals for the Sievers 500 RL^e On-Line TOC Analyzer from our Web site. Go to www.geinstruments.com, click Library, then Manuals, and make your selection.

Installation Instructions

Step 1 Unpack and inspect the Analyzer.

Step 2 Complete the identification records.

Step 3 Select a location for the Analyzer. Refer to the System Specifications chart.

The Sievers 500 RL^e On-Line TOC Analyzer is designed to be mounted on a wall or support stand. Allow a minimum of 5 cm clearance between the back of the Analyzer and the wall for heat dissipation; allow 30.5 cm of clearance on the sides, top, and bottom of the analyzer for the plumbing and electrical connections. Additionally, this clearance provides for the proper circulation for temperature and humidity control. When selecting the location, mount the Analyzer so that the display screen is approximately at eye-level.

The mounting hardware you provide must be able to support four times the weight of the unit. Thus, you should install mounting bolts capable of supporting 67.6 kg (148.8 lbs).

Avoid direct sunlight and extreme temperatures; operating at elevated temperatures greater than 40 °C (104 °F) prevents proper operation, and operating at low temperatures, 10 °C (50 °F), can cause errors in the measurements.

The Analyzer accessories kit contains two adhesive clips that can be used to store the double-ended screwdriver and John Guest fitting removal tool inside the Analyzer. Remove the paper backing from each clip and then firmly press the clip to the desired location inside the Analyzer. Do not attach the clips to any area that could affect analysis. Place the clips inside the door or on the right bulkhead.

Step 4 Install power and control wiring (Figure1).

Step 4 A Access the electrical enclosure.

Unlatch and open the door of the Analyzer. Locate and remove the plate over the electrical components, located on the far-left side. Remove this cover using a Phillip's (cross-head) screwdriver to loosen the two set screws.

Step 4 B Install a power supply.

This step should only be performed by a qualified electrician. The Analyzer requires an external source of AC power connected to the enclosure using a water-tight conduit connector. An external switch or circuit breaker should be installed near the Analyzer (clearly marked as the disconnecting device for the Analyzer).

Step 4 C Install the analog outputs and alarms.

The Analyzer offers four alarm and three 4-20 mA outputs, via three terminal blocks that are arranged horizontally in the electrical enclosure. Consult Tables 1 through 3 for a list of functions on each terminal block. The output and alarm connections should be installed by a qualified electrician.


Figure 1 - Wiring AC Power

English

Sievers 500 RL^e On-Line TOC Analyzers

Quick Start Guide

Warnings

See the separate "Warnings" document for additional important information.

The iOS port contains two sharp needles designed to pierce the septa of sample vials. Do not put fingers or inappropriate materials into the iOS port.

To avoid false TOC readings and possible damage to the Analyzer, always make sure the sample is flowing through the iOS System and the DI water reservoir is filled before starting analysis.

To avoid damaging the Analyzer and voiding the warranty, install the in-line filter on the sample inlet and replace the filter element as needed.

This symbol on the instrument indicates that the product does contain restricted substances included in China RoHS II. Please refer to the following table.

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Table of Hazardous Substances' Name and Concentration

部件名称 Component Name	有害物质 Hazardous Substances' Name					
	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr(VI))	多溴联苯 (PBB)	多溴二苯醚 (PBDE)
紫外灯 - UV Lamp	0	X	0	0	0	0
印刷电路板 - Printed circuit boards	X	0	0	0	0	0
本表格依据SJ/T11364的规范编制。 This table is prepared according to SJ/T 11364. 0: 表示该有害物质在该部件所有均质材料中的含量均在 GB/T 26572规定的限量要求以下 X: 表示该有害物质至少在该部件的某一均质材料中的含量超出 GB/T 26572规定的限量要求 • 此表所列数据为发布时所能获得的最佳信息 0: Indicates that hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in GB/T 26572. X: Indicates that t hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in GB/T 26572. • Data listed in the table represents best information available at the time of publication.						

Step 4 D Wiring the Ethernet cable.

Data from the Analyzer may be exported via Ethernet. Attach one end of the cable to the Ethernet port on the Analyzer's system board. Attach the other end of the cable to a data acquisition system. You must also activate the Modbus feature, as described in Chapter 4 of the 500 RL[®] O&M Manual.

Step 5 Install the printer, USB, and serial connections.

Connection ports for a printer cable, USB cable, and Serial cable are located on the top left side of the instrument. To view and access these ports, remove the cover plate by loosening the two thumb screws (Figure 2). Connect the appropriate cables now.


Figure 2 - Input / Output Connections

Step 6 Connect the sample inlet and outlet ports.

Before you begin, disable the the water flow to the instrument until the sample inlet system is completely installed and the Analyzer is ready to begin analysis.

1. Connect the 1/4" Teflon tubing with the in-line filter to the sample inlet on the iOS or Super iOS System or the Sample Inlet Block (depending on your Analyzer configuration). Tighten 1/4 turn past finger-tight with a 9/16" open-end wrench. Do not over-tighten the nut.
2. Connect the 3/4" OD waste line tubing to the waste outlet on the sample inlet system or the Sample Inlet Block (slide tubing over the barb fitting).
3. Place the hose clamp over the waste line and tighten to secure the connection to the waste outlet.
4. Route the waste tubing to an appropriate waste outlet. The waste is gravity drained. The waste tubing cannot be routed above the level of the waste outlet barb.
5. After water flow to the sample inlet system has been established, adjust the flow rate out of the waste line to between 50-300 mL/min. To adjust, turn the screw on the sample inlet system clockwise to decrease flow, and counter-clockwise to increase flow.

Step 7 Fill the DI water cartridge.

Fill the DI water cartridge with DI water, using the squeeze bottle provided in the accessories kit.

1. Locate the DI water cartridge at the center of the Analyzer (Figure 3).
2. Lift the spring clip (Figure 4) at the top of the DI water cartridge, and rotate the cartridge clockwise (as viewed from the top) so it slides out of the retaining hooks.
3. Remove the John Guest fittings labeled "B" and "D" from the barbs on the top of the DI cartridge.
4. Slide the water bottle nozzle into the inlet hole on the water cartridge.
5. Slide the DI water cartridge into the bracket and rotate it into the retaining clip until it snaps into place.
6. Put B fitting on D barb. Turn on the analyzer using the main power switch. This will start the DI pump. Run for 2 minutes.
7. Check for large air bubbles in all the tubing that connects to the DI water cartridge. If large air bubbles are present, tap or manipulate the tubing to work the air bubbles out. Put fitting B on barb B and D on D. Run for 1 minute.
8. Fill DI Water Cartridge to the line on the label (Figure 3).


Figure 3 - Location of DI Water

Figure 4 - Fill DI Water Reservoir

Step 8 Configure basic Analyzer settings

Before using the Analyzer, configure various basic settings. Refer to the 500 RL[®] On-Line O&M Manual for complete details.

- Step 8 A Enable password protection (optional).**
- Step 8 B Set the clock and time zone.**
- Step 8 C Name the Analyzer location (optional).**

Step 8 D Set the Analyzer mode.

Step 8 E Set up the data history.

Step 8 F Set up the printer (optional).

Step 8 G Export and print system settings.

Step 8 H Set Up Data I/O.

Step 9 Flush the Analyzer

Before placing the Analyzer into normal operation, rinse the Analyzer for 12 hours by running it in Rinse mode and then let it run in On-Line mode.

1. Select the Maintenance tab.
2. Press the Advanced button, and then press the Advanced Setup button.
3. Press the Auto Restart button. Make sure Rinse Down is set to On. If it is set to Off, press the Rinse Down button, and then press the On button.
4. Press the Back button, and then press the Rinse button.
5. The Analyzer will run in Rinse mode, and then will automatically switch into On-Line mode. Let the Analyzer run for 12 hours. No data will be reported for the Rinse mode.
6. After 12 hours, press the Stop Analysis button.
7. Open the Analyzer case and again visually inspect the Analyzer for leaks, especially around the DI water cartridge. If a leak is detected, make sure all fittings are tight and secure. Also confirm the water level in the DI water cartridge and, if necessary, turn off the Analyzer and add more water.

System Specifications*	
Linear range	0.02 ppb - 2,500 ppb TOC
Precision	≤1% RSD or 0.03 ppb TOC, whichever is greater
Accuracy	± 0.1 ppb or ± 5% of measurement, whichever is greater
Analysis mode	On-line; On-line Averaged; On-line Timed
Analysis time	Continuous mode: 6 minutes
Ozone Compatibility	50 ppb O ₃ continuous; 200 ppb O ₃ for 2 hours daily
Sample flow rate (nominal)	Analysis mode: 0.25 mL/min Fast Flush (between samples): 0.5 mL/min
External flow rate	Minimum 50 mL/min
Sample temperature	1 °C to 95 °C (34 °F to 203 °F)
Sample Pressure	Up to 100 psig
Interferences	Insensitive to organics heteroatoms
Calibration Stability	Typically stable for 12 months
Display Readout	3 significant digits
pH Range	5.5 - 8.0 in order to meet accuracy spec
Maximum Sample Conductivity	25 µS /cm at neutral pH
Instrument	
Power requirements	100-240 ±10% VAC, 100 watt, 50/60 Hz
Fuses	No user-replaceable fuses
Normal operating environment	Intended for indoor use only
Ambient temperature	10 °C to 40 °C (50 °F to 104 °F)
Maximum relative humidity	Up to 95%, noncondensing
Maximum altitude	3000 m (9,843 ft)
Inputs	One isolated binary input
Outputs	Serial (RS-232) ports, USB, parallel printer port, three 4-20 mA ports, four alarm ports, one Ethernet output
Installation/Overvoltage category	II (protects against transients present in Category II power)
Safety certificates	CE, ETL listed. Conforms to UL Std. 61010-1 Certified to CSA C22.2 No. 61010-1.
Pollution degree	2 (normally only non-conductive pollution)
Display	Backlit Quarter-VGA touchscreen display
Size	41.9 cm height x 48.3 cm width x 27.4 cm depth (16.5 in x 19.0 in x 10.8 in)
Weight	16.9 kg (37.2 lbs)
IP Rating	IP 45

* Stated analytical performance is achievable under controlled laboratory conditions that minimize operator and standards errors.

GE Power
Water & Process Technologies

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©2016. General Electric Company. Tutti i diritti riservati.